




# ***TDM ELISA***

# *General Information*

- TDM-ELISA consists of a colorimetric competitive enzyme immunoassay for the determination of HIV-1 protease inhibiting (PI) drugs in human plasma
- Correct Plasma concentrations of PI are needed for an efficient therapy and to avoid side effects
- PI's monitoring adds an important information in the management of HIV-positive patients

# *ELISA outline*

## *I. COMPETITION*




**Solid phase common to all TDM ELISA**

**Automation possible**

## *II. Detection Reaction*

**Chromogenic Substrate**



**Anti-drug antibody**

**Anti-rabbit antibody**

Hereafter is shown an outline of the preparation steps for a clear fresh plasma sample.

	<b>LPV</b>	<b>APV</b>	<b>NFV</b>	<b>RTV booster</b>
<b>Dilution in absolute methanol</b>	✓	✓	✓	
<b>Centrifugation</b>	✓		✓	
<b>Dilution in water</b>	✓		✓	
<b>Dilution in 30% methanol</b>	✓		✓	✓

Samples can be pre-centrifuged, aliquoted and stored at  $-20^{\circ}\text{C}$  avoiding repeated freezing and thawing

# Summary

Description and code	Sample Preparation	Assay procedure
<p>TDM-ELISA LPV</p> <p>Cod. 2678</p>	<p>100µl of plasma + 300µl of MetOH – Vortex for 10-15sec            Centrifugation : 10min at 10000 g in a microcentrifuge            Dilute 100µl of clear surnatant with 150µl of water.            Vortex            Add to 100µl of the diluted surnatant 150µl of MetOH 30%            Vortex            Use 20 µl for testing</p>	<p>20µl of pre-treated sample            +            80µl of specific enzyme            +            100 µl specific Antiserum</p> <p><u>Incubate 1h at RT</u></p> <p>Wash wells 5 times with Washing Buffer</p> <p>Add 200µl of chromogenic solution ( TMB) diluted in the development solution</p> <p><u>Incubate 30min (RT in the dark)</u></p> <p>Add 50µl Of stop solution</p> <p>Read at 450nm on a microplate reader</p>
<p>TDM-ELISA RTV booster</p> <p>Cod. 7678</p>	<p>10µl of + 990 µl of MetOH 30%</p> <p>Vortex</p> <p>Use 20 µl for testing</p>	
<p>TDM-ELISA APV / Fos APV</p> <p>Cod. 1678</p>	<p>10µl of plasma + 990 µl of sample diluent</p> <p>Vortex</p> <p>Use 20 µl for testing</p>	
<p>TDM-ELISA NFV</p> <p>Cod. 3678</p>	<p>100µl of plasma + 300µl of MetOH – Vortex for 10-15sec            Centrifuge for 10min at 10000 g in a microcentrifuge            Dilute 100 µl of clear surnatant with 50 µl of distilled water.            Vortex            Add to 10 µl of the diluted surnatant 990µl of 30%.Methanol            Vortex            Use 20 µl for testing</p>	

# *Applications*

**This technology allows the determination of PI concentration in plasma samples of patients in HAART treatment**

**TDM-ELISA is especially useful in the following conditions:**

- ★ **Pregnant women and children**
- ★ **Patients with suspected significant pharmacologic interactions which might induce dose-related adverse effects or which could cause a reduction of treatment efficacy**
- ★ **Physiopathological conditions significantly modifying liver or kidney's function and inducing alterations in absorption, distribution, metabolism or clearance of the drug**
- ★ **Patients with viral strains resistant to treatment and needing different dosage schedules whose efficacy has not yet been proved in clinical trials**
- ★ **Dose-related toxicity**
- ★ **Lack of virologic response in “naive” patients**

# **IMPORTANCE OF TDM - ELISA**

**A normal TDM result is always an important information for the clinician because it allows to exclude abnormal drug levels as cause of toxic symptoms or uneffective therapy**

# Data from Medizinische Laboratorien Düsseldorf (Germania) - dr Kuschak

## TDM-ELISA LPV

Table 1 - Summary

7 control samples were used from Ring Test at known concentration for 7-10 times on the same plate

# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
8240	7649	7516	8842	8681	8681	7709	8681	6962		8195	8106,8	668,5	8,2	-1,1
# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
7085	7827	7440	7287	7709	8053	6387	7515	7929		7990	7470,2	511,5	6,8	-6,5
# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
3155	3171	2863	3194	3185	3521	3528	3478	3444	3419	2930	3295,8	216,1	6,6	12,5
# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
1569	1808	1831	1822	1776	1873	1822	1799	1849	1901	1639	1805,0	90,4	5,0	10,1
# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
785	792	864	1011	841	854	892	869	900	864	737	867,2	63,0	7,3	17,7
# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
3057	3329	3049	3027	3393	3777	3768	3486	3328	3695	3296	3390,9	290,9	8,6	2,9
# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
15950	16080	15860	16210	16760	17760	18210				16214	16690	939,8	5,6	2,9

# TDM-ELISA APV

## Summary table 2

Evaluation on 4 Ring Test samples repeated 9 - 10 times (intra-assay)

# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
5131	7394		6484	6862	7510	7451	5756	6183	3376	5720	6238,6	1347,6	21,6	9,07
# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
3887	2489	2409	3077	4357	3560		2660	3609	2896	3600	3216,0	674,1	21,0	-10,67
# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
3039	1984	3361	2930	2711	3248	2700	2516	3659	2023	2860	2817,1	547,0	19,4	-1,5
# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
1596	2010	1505	2010	1721	1991	2120	2581	2085	2023	1800	1964,2	304,3	15,5	9,1

# TDM-ELISA RTV Booster

Table 3 - Summary

7 Ring Test samples at known concentration were use for 10 times on the same plate (intra-assay)

# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
1039	1127	1250	927	986	1067	1144	1177	1144	1160	1100	1102	100,4	9,1	0,2
# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
492	447	511	477	500	614	492	573	532	532	550	517,0	51,0	9,9	-6,0
# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
197	187	186	218	169	187	209	210	190	210	220	196,3	15,2	7,7	-10,8
# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
1025	1122	1210	1115	1115	1161	971	1087	947	965	1172	1071,8	90,0	8,4	-8,5
# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
77,4	88	95,7	85,9	85,3	84,8	49,7	91,9	77,7	112	110	84,8	15,8	18,7	-22,9
# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
55,2	51,8	48,5	46,1	83	48,1	50,2	50,6	51,3	50,6	73,3	53,5	10,6	19,9	-27,0
# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9	# 10	Target Value	Mean	St Dev	CV	Inaccuracy
30,5	37,7	45,9	39,7	38,5	39,4	44,3	43,3	43,6	40,8	52,4	40,3	4,7	11,6	-23,0

# *Results of the evaluation*

## *Accuracy and Precision*

The data obtained with TDM-ELISA are equivalent to the data which characterize the performance of HPLC

( Drose J.A. et. Alt., 32, 287-291 – JAIDS 2003)

# *TDM-ELISA*

## *Specificity*

- ★ No cross-reactivity has been observed of each TDM with respect to the other peptidic PI  
(APV, LPV, RTV, NFV e M8, ATV, IDV, SQV)  
and with **TIPRANAVIR** (new drug of a new class of non-peptidic PI)
- ★ No cross-reactivity has been observed with other drugs used in HAART therapy  
(NNRTIs: Nevirapine, Efavirenz; NRTIs: Zidovudine, Abacavir and Lamivudine)

## *Detection limit in matrix*




TDM - ELISA	ng/ml
LPV	123
APV	400
RTV	40

## *Measuring Range*

- ★ In agreement with the American Guidelines (US Therapy – DHHS ) of October 2006 the kits have measuring ranges based on the significant therapeutic limit

## Dynamic measuring ranges of the different kits

Drug	Quantification Range Accuracy within 20 % (concentrations: $\mu\text{g} / \text{ml}$ )
LPV	1 – 8
NFV	0,8 - 3
APV	0,4 - 5
RTV booster	0,09 – 0,9
RTV	0,5 - 5



# Conclusions

---

ELISA TDM for their low cost and simplicity can be a reliable and useful method for routine monitoring of TDM in any laboratory

## *Availability*

### ★ Available Products :

<b>TDM-ELISA APV</b>	<b>cod. 1678</b>
<b>TDM-ELISA LPV</b>	<b>cod. 2678</b>
<b>TDM-ELISA NFV</b>	<b>cod. 3678</b>
<b>TDM-ELISA RTV booster</b>	<b>cod. 7678</b>
<b>TDM-ELISA ATV*</b>	<b>cod. 5678</b>